

Last update: 25 October 2023

REACH Regulation (EC) 1907/2006

Article 33 statement

To professional/commercial users or distributors of Schiller devices

One or more articles shipped by SCHILLER AG contain substances which were included in the European Chemicals Agency's list of *Substances of Very High Concern* (SVHC), available at <https://www.echa.europa.eu/candidate-list-table>.

According to *Regulation (EC) 1907/2006 Article 33* ("REACH"), suppliers of articles containing a substance in the SVHC list in a concentration higher than 0.1 % by weight are required to communicate information regarding safe use of the substances. Please note that inclusion of a substance in the SVHC list does not restrict its placement on the market or its use.

SCHILLER AG will continue to monitor the regulatory requirements related to chemical substances and will update this document should the need arise.

The list of substances and articles concerned can be found in the annex of the present document.



SCHILLER

The Art of Diagnostics

Annex: List of substances in SVHC list used in Schiller devices

The substances listed below may be contained in parts of the specified devices in a concentration higher than 0.1% by weight. Devices and parts that do not contain SVHCs are not listed. No specific instructions on safe use are required as the substances in question will not be released in significant amounts under normal or reasonably foreseeable conditions of use.

Substance Name	Dodecachloropentacyclo[12.2.1.16,9.02,13.05,10]octadeca-7,15-diene or Dodecachlorodimethanodibenzocyclooctane (CAS 13560-89-9)	
Use	Flame-retardant additive in cable jackets	
Article number and name of devices containing the substance	0.775000	ARGUS PRO LifeCare 2
	0A.107000	CARDIOVIT AT-1 G2
	0.080000	CARDIOVIT AT-101
	0.070000	CARDIOVIT AT-102
	0.075000	CARDIOVIT AT-102 plus
	0A.108000	CARDIOVIT AT-102 G2
	0.030000	CARDIOVIT CS-200
	0.030100	CARDIOVIT CS-200 ErgoSpiro
	0.035000	CARDIOVIT CS-200 Excellence
	0A.102000	CARDIOVIT CS-200 Touch Box
	0A.800000	Diagnostic Station DS20
	0.760215	ARGUS PB-1000LP Merge 4xIBP
	0.760400	ARGUS PB-1000LP Cath Med 4xIBP



SCHILLER

The Art of Diagnostics

Substance Name	Octamethylcyclotetrasiloxane (CAS 556-67-2)	
Use	Manufacturing process related silicone component used in splitter cable and Y-cable box	
Article number and name of devices containing the substance	0.790000	ARGUS LCX
	0.760215	ARGUS PB-1000LP Merge 4xIBP
	0.760255	ARGUS PB-1000LP 4xIBP
	0.760400	ARGUS PB-1000LP Cath Med 4xIBP
	0.770000, 0.770200, 0.771000	ARGUS PB-2000LP
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000

Substance Name	4,4'-isopropylidenediphenol (CAS 80-05-7)	
Use	Precursor to plastics used in printed circuit boards	
Article number and name of devices containing the substance	0A.107000	CARDIOVIT AT-1 G2
	0A.110000	CARDIOVIT AT-180
	0.035000	CARDIOVIT CS-200 Excellence
	0.035100	CARDIOVIT CS-200 Excellence ErgoSpiro
	0A.900000	FRED easyport plus
	0A.500000	SPIROVIT SP-1 G2
	1A.702100	Tempus LS
	1A.702200	Tempus LS MANUAL



SCHILLER

The Art of Diagnostics

Substance Name	Lead (CAS 7439-92-1)	
Use	Copper alloy component of connector pins or EEE components	
Article number and name of devices containing the substance	0.780000, 0.781000	ARGUS LCM plus/basic
	0.775000	ARGUS PRO LifeCare 2
	0A.107000	CARDIOVIT AT-1 G2
	0.080000	CARDIOVIT AT-101
	0.070000	CARDIOVIT AT-102
	0A.108000	CARDIOVIT AT-102 G2
	0A.110000	CARDIOVIT AT-180
	0.310000	BP-200 plus
	0.290000	BR-102 plus
	0A.111000	CARDIOVIT CS-104
	0.035000	CARDIOVIT CS-200 Excellence
	0.035100	CARDIOVIT CS-200 Excellence ErgoSpiro
	0A.103000	CARDIOVIT CS-200 Office Blue
	0A.103010	CARDIOVIT CS-200 Office Blue – PAKET
	0A.105000	CARDIOVIT PRO PC System
	0A.102000	CARDIOVIT CS-200 Touch Box
	0A.107100	CT8000i-2
	0A.108100	CT8000P-2
	0A.800000	Diagnostic Station DS20
	0A.400000	EASY PULSE
	0.900000	FRED easyport
	0A.900000	FRED easyport plus
	0A.106000	CARDIOVIT FT-1
	0A.112000	CARDIOVIT FT-2
	1A.306000	medilogAR Holter
	0A.100000	CARDIOVIT MS-2015
	0.760000	ARGUS PB-1000LP
	0.760215	ARGUS PB-1000LP Merge 4xIBP
0.760255	ARGUS PB-1000LP 4xIBP	
0.760400	ARGUS PB-1000LP Cath Med 4xIBP	



SCHILLER

The Art of Diagnostics

	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000
	0A.500000	SPIROVIT SP-1 G2
	1A.702100	Tempus LS
	1A.702200	Tempus LS MANUAL
	1A.901000	Trainer FRED easyport plus

Substance Name	Imidazolidine-2-thione (CAS 96-45-7)	
Use	Vulcanization accelerator in NIBP pump O-rings	
Article number and name of devices containing the substance	0.775000	ARGUS PRO LifeCare 2
	0.780000, 0.781000	ARGUS LCM plus/basic
	0.310000	BP-200 plus
	0A.800000	Diagnostic Station DS20
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000
	0.311000	QT-1
	3.930430	Welch Allyn 1500 Patient Monitor



SCHILLER

The Art of Diagnostics

Substance Name	2-methyl-1-(4-methylthiophenyl)-2-morpholinopropan-1-one (CAS 71868-10-5)	
Use	EEE components	
Article number and name of devices containing the substance	0A.110000	CARDIOVIT AT-180
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000



SCHILLER

The Art of Diagnostics

Substance Name	Azodicarbonamide (CAS 123-77-3)	
Use	Blowing agent in adhesive tape	
Article number and name of devices containing the substance	0.780000, 0.781000	ARGUS LCM plus/basic
	0.775000	ARGUS PRO LifeCare 2
	0.170000	CARDIOVIT AT-10
	0.080000	CARDIOVIT AT-101
	0.290000	BR-102 plus
	0.083000	CARDIOVIT AT-101 easy
	0.084000	CARDIOVIT AT-101 Tele
	0A.800000	Diagnostic Station DS20
	0.090000	CARDIOVIT MS-2010
	1.300000, 1.300001, 1.300010, 1.300011	MICROVIT MT-101
	0.091000	MyCardioPad
	0.081000	P8000
	0.085000	P8000 Tele
0.540000	SPIROVIT SP-1	

Substance Name	Lead titanium trioxide (CAS 12060-00-3)	
Use	EEE components	
Article number and name of devices containing the substance	0A.110000	CARDIOVIT AT-180
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000



SCHILLER

The Art of Diagnostics

Substance Name	Cadmium hydroxide (CAS 21041-95-2)	
Use	Use in battery pack	
Article number and name of devices containing the substance	0.040000, 0.040100, 0.040200	CARDIOVIT AT-104

Substance Name	Bis(2-(2-methoxyethoxy)ethyl) ether (CAS 143-24-8)	
Use	EEE components	
Article number and name of devices containing the substance	0A.900000	FRED easyport plus

Substance Name	1-methyl-2-pyrrolidone (CAS 872-50-4)	
Use	EEE components	
Article number and name of devices containing the substance	0A.900000	FRED easyport plus

Substance Name	2,2',6,6'-tetrabromo-4,4'-isopropylidenediphenol (CAS 79-94-7)	
Use	EEE components	
Article number and name of devices containing the substance	0A.110000	CARDIOVIT AT-180
	0A.107000	CARDIOVIT AT-1 G2
	0A.108000	CARDIOVIT AT-102 G2
	0.290000	BR-102 plus
	0A.106000	CARDIOVIT FT-1
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000
	0A.103100, 0A.103110, 0A.103200	CARDIOVIT CS-200 Office
	0A.109000, 0A.111000, 0A.111010	CARDIOVIT CS-104
	0.900000	FRED easyport
	0A.900000	FRED easyport plus



SCHILLER

The Art of Diagnostics

	1A.305000	MS-12 USB
	0A.500000	SPIROVIT SP-1 G2

Substance Name	Diboron trioxide (CAS 1303-86-2)	
Use	EEE components	
Article number and name of devices containing the substance	0.030000	CARDIOVIT CS-200
	0.040000, 0.040100, 0.040200	CARDIOVIT AT-104
	0.060000	CARDIOVIT AT-110
	0.170000	CARDIOVIT AT-10
	1A.701007, 1A.701010, 1A.701101, 1A.701105, 1A.701107	ARGUS PB-3000
	0A.109000, 0A.111000, 0A.111010	CARDIOVIT CS-104
	0A.103100, 0A.103110, 0A.103200	CARDIOVIT CS-200 Office
	0.030100	CARDIOVIT CS-200 ErgoSpiro

Substance Name	Dodecamethylcyclohexasiloxane (CAS 540-97-6)	
Use	EEE components	
Article number and name of devices containing the substance	2.000202	Kalibrations-Kit Hans Rudolph