

BP-200 plus

Automated blood pressure monitor for cardiac exercise testing



SCHILLER
The Art of Diagnostics

BP-200 plus

- ∴ Seamless integration: automated communication with your SCHILLER exercise ECG device reduces the risk of transmission errors.
- ∴ Non-exercise BP mode: allows BP measurements during patient setup and recovery without ECG signal.
- ∴ Verified measurements: Korotkoff sounds are displayed.



BP AND SPO₂ MEASUREMENTS FOR THE EXERCISE TEST ENVIRONMENT

The BP-200 plus blood pressure monitor is designed to give you the confidence to know that mechanical vibration, patient movement and observer variability does not interfere with an accurate reading. Your focus can be where it matters most – on your patient. The hands-free interface and SpO₂ measurement option make the BP-200 plus a comprehensive device for diagnostic test centres. The BP-200 plus is also ideal in research environments because raw data can be saved on an SD card.

SCHILLER K-SOUND ANALYSIS FOR ACCURATE MEASUREMENTS

The BP-200 plus excels where other BP monitors fail, by using SCHILLER K-Sound Analysis (SKA). SKA was developed as a result of SCHILLER's pioneering work in 24-hour ambulatory blood pressure monitoring. It was this experience that provided SCHILLER's research and development team with invaluable insights into the relationship between many of the variables involved in determining the blood pressure of a patient in motion. These variables include, but are not limited to, K-sounds, overall body shape, foot strike characteristics, arm movements, along with treadmill and ergometer noise. With SCHILLER K-Sound Analysis, Korotkoff sounds are distinguished from artefact noise so that blood pressures are measured accurately with greatly reduced false measurements.

AUTOMATED COMMUNICATION WITH SCHILLER EXERCISE ECG SYSTEM

Recording an accurate and reliable BP measurement or SpO₂ value at the correct intervals during an exercise ECG test can be difficult and stressful for clinicians and lab technicians. To increase measurement reliability, the BP-200 plus adopts the programmed intervals from the SCHILLER exercise ECG device, ensuring that measurements are taken at the same point during each stage. For additional measurement accuracy and stability, the QRS trigger from the exercise ECG system may be used.

Systolic BP, diastolic BP, mean arterial pressure, SpO₂ (option), peripheral pulse and double product (sys × HR)

5 buttons for easy programming and operation



Display for Korotkoff sounds or trigger waves

SCHILLER BP cuff is designed for best comfort during exercise ECGs





Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
 sales@schilleramericas.com
 www.schilleramericas.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
 sales@schilleramericas.com
 www.schilleramericas.com



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 sales@schiller.com.my
 www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99 0
 Fax +43 732 757 000
 sales@schiller.at
 www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co., Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
 info@schillerchina.com
 www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
 info@schillerzg.hr
 www.schiller.hr



France
SCHILLER Médical
 F-67160 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
 info@schiller.fr
 www.schiller.fr



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
 contact@schillerfrance.fr
 www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
 info@schillermed.de
 www.schillermed.de



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai – 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
 sales@schillerindia.com
 www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
 schiller@schiller.pl
 www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
 119049 Moscow, Russia
 Phone +7 (495) 970 11 33
 mail@schiller.ru
 www.schiller.ru



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
 info@schiller.rs
 www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
 info@schiller.si
 www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
 schiller@schiller.es
 www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8912 Obfelden
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
 sales@schiller-reomed.ch
 www.schiller-reomed.ch



Turkey
SCHILLER TÜRKIYE
 Okmeydani-Sisli – Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
 info@schiller.com.tr
 www.schiller-turkiye.com



United Kingdom
SCHILLER UK Ltd.
 Anstruther, Fife KY10 3H
 Phone +44 16 190 72 565
 sales@schilleruk.com
 www.schilleruk.com

Device availability in your market is subject to regulatory approval.

Article no.: 2.500248
 2021/12, EN | Subject to changes.



SCHILLER
 The Art of Diagnostics