

CARDIOVIT AT-102 G2

Easy, smart and state-of-the-art



SCHILLER
The Art of Diagnostics

CARDIOVIT AT-102 G2

The new generation ECG: CARDIOVIT AT-102 G2

- ✦ Best signal quality and diagnostic algorithms for adult and paediatric patients
- ✦ Fast bidirectional Wi-Fi communication with HIS/EMR
- ✦ Spirometry and basic exercise ECG optionally available



SpiroScout SP plus, for spirometry tests with the CARDIOVIT AT-102 G2

STATE-OF-THE-ART TECHNOLOGY AND CONNECTIVITY

The CARDIOVIT AT-102 G2 features an 8" colour display with touch function keys. Its sealed alphanumeric keyboard ensures easy cleaning, making it highly suitable for the daily clinical routine. Smart battery technology guarantees more than 8 hours of ECG recording.

Bidirectional communication allows for easy data access and fast transmission of all reports to EMR/HIS systems in PDF, DICOM or HL7 format. Standard Wi-Fi and LAN connectivity, combined with ECG review on the display and a large memory, supports paperless workflows and cost saving. Customizable printouts are available on the internal high-resolution A4 thermal printer.

HIGH EFFICIENCY IN ANY WORKFLOW

ECG acquisition has never been easier and faster: Enter or download patient data with PDQ or worklist, acquire ECG, save, print and transmit data within twenty seconds using only one navigation key.




EASY 1-2-3 STEPS

- 1: Enter or download patient data
- 2: Check signal quality and perform ECG
- 3: Review ECG before storing, exporting or printing

THE CARDIOVIT AT-102 G2 DIAGNOSTIC TOOLS

- ✦ Resting ECG with the SCHILLER ECG analysis program ETM as well as ETM Sport, the automatic interpretation for athletes.
- ✦ Resting rhythm recording up to 20 minutes.
- ✦ Optimal signal quality for paediatric ECGs.
- ✦ Basic exercise ECG (optional) with a large selection of protocols for treadmills and ergometers.
- ✦ SCHILLER's Culprit Coronary Artery Algorithm (CCAA, optional) locates the obstruction site in the coronary artery and thereby the size of the cardiac area at risk.



The CARDIOVIT AT-102 G2 is also available with a trolley.

ADD SPIROMETRY TO YOUR CARDIOVIT AT-102 G2

With the SpiroScout SP plus, the CARDIOVIT AT-102 G2 also becomes a spirometer. Thanks to GANSHORN's ultrasound technology, spirometry recordings are highly accurate and the sensor is maintenance- and calibration-free. The disposable ScoutTubes ensure maximum hygiene and patient safety. ScoutTubes are CO₂ neutral.

The spirometry option offers FVC, SVC and MVV with a large selection of reference standards, pre and post tests, as well as interpretation of the FVC results.

The device provides feedback on the procedure's quality during the test and after each trial.

Complete spirometry results can be printed on thermal paper as well as exported to EMR/HIS systems.



An incentive screen is included for trials with children.



Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com



France
SCHILLER Médical
F-67160 Wissembourg
Phone +33 3 88 63 36 00
info@schiller.fr
www.schiller-medical.fr



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
contact@schillerfrance.fr
www.schiller-france.com



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
info@schiller.si
www.schiller.si



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
sales@schiller.com.my
www.schiller-asia.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
info@schillermed.de
www.schillermed.de



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
schiller@schiller.es
www.schiller.es



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 99 0
sales@schiller.at
www.schiller.at



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 059, India
Phone +91 22 6152 3333 / 2920 9141
sales@schillerindia.com
www.schillerindia.com



Switzerland
SCHILLER Schweiz AG
CH-8912 Obfelden
Phone +41 44 744 30 00
sales@schiller-schweiz.ch
www.schiller-schweiz.ch



China
Alfred Schiller (Beijing) Medical Technology
Co., Ltd.
100102 Beijing, China
Phone +86 10 84565453
info@schillerchina.com
www.schillermedical.cn



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
schiller@schiller.pl
www.schiller.pl



Turkey
SCHILLER TÜRKİYE
Okmeydani-Sisli - Istanbul
Phone +90 212 210 8681
info@schiller.com.tr
www.schiller-turkiye.com



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
info@schiller.hr
www.schiller.hr



Russia & C.I.S.
AO SCHILLER.RU
119049 Moscow, Russia
Phone +7 (495) 970 11 33
mail@schiller.ru
www.schiller.ru



United Kingdom
SCHILLER UK Ltd.
Bromborough, Wirral,
Merseyside, CH62 3RJ
Phone +44 1333 312150
sales@schilleruk.com
service@schilleruk.com
www.schilleruk.com

SCHILLER reserves the right to make changes in design, specifications and features shown or discontinue the product described at any time without notice or obligation. Images may be representative. Some products may not be regulatory cleared or released for sale in all markets. Contact your SCHILLER representative for the most current information and availability.

