


CARDIOVIT CS-104

The exercise ECG of the future in an integrated system




SCHILLER
The Art of Diagnostics

CARDIOVIT CS-104

Everything you need in
one system:

- ∴ Wireless ECG acquisition
- ∴ Optional ultrasound spirometry with SpiroScout SP plus
- ∴ Integrated workstation for data review



ALL POSSIBILITIES IN A SINGLE SYSTEM

- ∴ Resting ECG with the SCHILLER ECG analysis program ETM as well as ETM Sport, the automatic interpretation for athletes
- ∴ Resting rhythm recording up to 20 minutes, with optional arrhythmia detection
- ∴ ECG Framer: 10-second resting ECG from resting rhythm
- ∴ Exercise ECG with optional arrhythmia detection
- ∴ Vector ECG (optional)
- ∴ Spirometry, based on the ultrasound technology of the SpiroScout SP plus (optional)



WORKSTATION FOR DATA REVIEW

The CARDIOVIT CS-104 provides extensive review features for both data locally acquired and recordings from all SCHILLER devices connected to the SCHILLER Server. Bi-directional connectivity (including worklist and PDQ) to EMR, HIS and PACS systems is given.

WIRELESS ECG

The wireless module MS-12 blue greatly simplifies ECG acquisition, allowing the user to freely operate in any environment and minimizing motion artefacts. Whether laying on the bed or running on the treadmill, the patient has the greatest possibility to move and will hardly notice this small and light device (115 g).

SWISS QUALITY - READY-TO-USE

With the CARDIOVIT CS-104, SCHILLER offers a complete ECG system, ready for use in practices, clinics and hospitals.

The unit includes a PC, a monitor, a trolley and a wireless (MS-12 blue) or wired (MS-12 USB) acquisition module. Optionally, the CARDIOVIT FT-1 can be used for ECG acquisition. Everything is pre-installed and tested by SCHILLER, ready to be used in the clinical practice.

CARDIOVIT CS-104 is also available as resting ECG PC software and spirometry PC software.

Ultrasound spirometer
SpiroScout SP plus

Intuitive and customizable user interface

Soft handle positioned at an adjustable height

Small but powerful wireless ECG acquisition module with on-screen signal view and hook-up advisor

Transparent bin to store and transport medical items

Easy-glide wheels to smoothly maneuver the trolley





Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com



France
SCHILLER Médical
 F-67160 Wissembourg
 Phone +33 3 88 63 36 00
info@schiller.fr
www.schiller.fr



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
info@schiller.rs
www.schiller.rs



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
contact@schillerfrance.fr
www.schiller-france.com



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
info@schiller.si
www.schiller.si



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
sales@schiller.com.my
www.schiller-asia.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
info@schillermed.de
www.schillermed.de



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
schiller@schiller.es
www.schiller.es



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99 0
sales@schiller.at
www.schiller.at



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai – 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
sales@schillerindia.com
www.schillerindia.com



Switzerland
SCHILLER-Reomed AG
 CH-8912 Obfelden
 Phone +41 44 744 30 00
sales@schiller-reomed.ch
www.schiller-reomed.ch



China
Alfred Schiller (Beijing) Medical Equipment Co., Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
info@schillerchina.com
www.schillermedical.cn



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
schiller@schiller.pl
www.schiller.pl



Turkey
SCHILLER TÜRKIYE
 Okmeydani-Sisli – Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
info@schillerzg.hr
www.schiller.hr



Russia & C.I.S.
AO SCHILLER.RU
 119049 Moscow, Russia
 Phone +7 (495) 970 11 33
mail@schiller.ru
www.schiller.ru



United Kingdom
SCHILLER UK Ltd.
 Anstruther, Fife KY10 3H
 Phone +44 16 190 72 565
sales@schilleruk.com
www.schilleruk.com

Device availability in your market is subject to regulatory approval.

