



CARDIOVIT CS-104

The integrated stress test system




SCHILLER
The Art of Diagnostics

CARDIOVIT CS-104

Everything you need in
one system:

- ❖ Wireless ECG acquisition
- ❖ Optional ultrasound spirometry with SpiroScout SP plus
- ❖ Integrated workstation for data review



Exercise stress ECG with CARDIOVIT FT-1

ALL POSSIBILITIES IN A SINGLE SYSTEM

- ❖ Exercise stress test with optional arrhythmia detection
- ❖ Pharmacological stress test
- ❖ Resting ECG with the SCHILLER ECG analysis program ETM as well as ETM Sport, the automatic interpretation for athletes
- ❖ Resting rhythm recording up to 60 minutes, with optional arrhythmia detection
- ❖ ECG Framer: 10-second resting ECG derived from resting rhythm

- ❖ Vector ECG (optional)
- ❖ Spirometry, based on the ultrasound technology of the SpiroScout SP plus (optional)

WORKSTATION FOR DATA REVIEW

The CARDIOVIT CS-104 provides extensive review and editing features as well as tools for final report preparation. Both locally acquired data and recordings from other SCHILLER resting ECG devices and exercise stress test systems connected to the SCHILLER Server can be examined. Bi-directional connectivity (including worklist and PDQ) to EMR, HIS, and PACS systems is given.

WIRELESS ECG

The optional wireless module MS-12 BT greatly simplifies ECG acquisition, allowing the user to operate freely in any environment and minimising motion artefacts. Whether biking or running on a treadmill, the patient can move freely and will barely notice this small, light-weight device (125 g).

SWISS QUALITY – READY-TO-USE

The system includes a PC, a monitor, a trolley, and a wireless (MS-12 BT) or wired (MS-12 USB) acquisition module. Optionally, CARDIOVIT FT-1 can be used for ECG acquisition. Everything is pre-installed and tested by SCHILLER, ready to be used in the clinical practice.

CARDIOVIT CS-104 is also available as PC software for resting ECG, with optional stress test and spirometry add-ons.



Ultrasound spirometer
SpiroScout SP plus

BP-200 plus (optional)
for blood pressure and
SpO₂ measurement
during the stress test

Intuitive and customi-
sable user interface:

- ❖ Choose from different ECG views
- ❖ Create your preferred workflows
- ❖ Define search filters and report templates

Choose between three
acquisition modules:

- ❖ MS-12 BT: Small but powerful wireless ECG acquisition module with on-screen signal view and hook-up advisor
- ❖ MS-12 USB: Easy-to-use plug-and-play ECG acquisition module
- ❖ CARDIOVIT FT-1: This SCHILLER standalone resting ECG device may also be used as an acquisition module

Working height
adjustable from 88 to
122 cm with a spring
mechanism

Add spirometry to
your stress test
system, with the
SpiroScout SP plus
(optional)

Transparent bin to
store and transport
medical items

Easy-glide wheels to
smoothly manoeuvre
the trolley





Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 99 0
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Technology Co., Ltd.
100102 Beijing, China
Phone +86 10 84565453
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
info@schiller.hr
www.schiller.hr



France
SCHILLER Médical
F-67160 Wissembourg
Phone +33 3 88 63 36 00
info@schiller.fr
www.schiller-medical.fr



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
info@schillermed.de
www.schillermed.de



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 059, India
Phone +91 22 6152 3333 / 2920 9141
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
119049 Moscow, Russia
Phone +7 (495) 970 11 33
mail@schiller.ru
www.schiller.ru



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER Schweiz AG
CH-8912 Obfelden
Phone +41 44 744 30 00
sales@schiller-schweiz.ch
www.schiller-schweiz.ch



Turkey
SCHILLER TÜRKİYE
Okmeydanı-Sisli - Istanbul
Phone +90 212 210 8681
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK Ltd.
Bromborough, Wirral,
Merseyside, CH62 3RJ
Phone +44 1333 312150
sales@schilleruk.com
service@schilleruk.com
www.schilleruk.com

SCHILLER reserves the right to make changes in design, specifications and features shown or discontinue the product described at any time without notice or obligation. Images may be representative. Some products may not be regulatory cleared or released for sale in all markets. Contact your SCHILLER representative for the most current information and availability.

