


CARDIOVIT CS-104 SPIRO


Networked PC spirometry with ultrasonic technology



CARDIOVIT CS-104 SPIRO

CARDIOVIT CS-104 Spiro is the all-in-one solution in lung function testing. A desktop spirometry solution for large and small healthcare providers:

- Offers networking options and SEMA integration
- Ultrasound technology: calibration-free, no warm-up time, no moving parts



Incentive screen for animation during spirometry

NETWORKED SPIROMETRY

CARDIOVIT CS-104 Spiro provides the advantages of network connectivity:

- Connected as required via the connectivity solution SCHILLER Server.
- Use the SCHILLER data management system SEMA. Results of pulmonary function tests can be transmitted and stored in SEMA for review and editing.

ULTRASOUND TECHNOLOGY

Thanks to GANSHORN's ultrasound technology, recordings are highly accurate, and the sensor is maintenance- and calibration-free. The disposable Scout-Tubes ensure maximum hygiene and patient safety.

REAL-TIME INFORMATION

The software helps to ensure that the patient achieves a good result and that the measurement is valid: information is displayed in real time, feedback on the procedure's quality is provided after each trial, and an incentive screen is available for trials with children.

WIDE RANGE OF MEASUREMENTS

The spirometer offers FVC, SVC, and MVV. Moreover, Pre and Post tests are possible. The CARDIOVIT CS-104 Spiro automatically interprets the results and provides information to support diagnosis.

PERFORM SPIROMETRY WHEREVER YOU ARE

The CARDIOVIT CS-104 Spiro offers the major advantage of performing spirometry wherever you are, with a portable PC.

DISPOSABLE AND SUSTAINABLE

To ensure a sustainable product, SCHILLER makes a contribution to climate protection for the ScoutTubes and the SpiroDef ScoutTube filters.


Engaged for Impact
Product
myclimate.org/01-20-162037



Connected to your infrastructure via SCHILLER Server

Compliance with the latest ATS/ERS Guidelines

Automatic interpretation of the results


Review and edit data with the SCHILLER data management system SEMA

FVC, SVC, MVV, Pre/Post testing, wide range of predicted normal values

Providing feedback on the procedure's quality

Disposable ScoutTubes for maximum hygiene and patient safety

Ultrasonic sensor: calibration-free and maintenance-free





Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schilleramericas.com
www.schilleramericas.com



France
SCHILLER Médical
F-67160 Wissembourg
Phone +33 3 88 63 36 00
info@schiller.fr
www.schiller.fr



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
contact@schillerfrance.fr
www.schiller-france.com



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
info@schiller.si
www.schiller.si



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
sales@schiller.com.my
www.schiller-asia.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
info@schillermed.de
www.schillermed.de



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
schiller@schiller.es
www.schiller.es



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 99 0
sales@schiller.at
www.schiller.at



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai – 400 059, India
Phone +91 22 6152 3333 / 2920 9141
sales@schillerindia.com
www.schillerindia.com



Switzerland
SCHILLER Schweiz AG
CH-8912 Obfelden
Phone +41 44 744 30 00
sales@schiller-schweiz.ch
www.schiller-schweiz.ch



China
Alfred Schiller (Beijing) Medical Technology
Co., Ltd.
100102 Beijing, China
Phone +86 10 84565453
info@schillerchina.com
www.schillermedical.cn



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
schiller@schiller.pl
www.schiller.pl



Turkey
SCHILLER TÜRKIYE
Okmeydani-Sisli – Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
info@schiller.hr
www.schiller.hr



Russia & C.I.S.
AO SCHILLER.RU
119049 Moscow, Russia
Phone +7 (495) 970 11 33
mail@schiller.ru
www.schiller.ru



United Kingdom
SCHILLER UK Ltd.
Anstruther, Fife KY10 3H
Phone +44 1333 312150
sales@schilleruk.com
www.schilleruk.com

SCHILLER reserves the right to make changes in design, specifications and features shown or discontinue the product described at any time without notice or obligation. Images may be representative. Some products may not be regulatory cleared or released for sale in all markets. Contact your SCHILLER representative for the most current information and availability.

