


# FRED® EASYPORT®

The world's first pocket defibrillator



# FRED® EASYPORT®

- ∴ Extremely light:  
only 490 grams (incl. battery)
- ∴ Extremely small:  
133 × 126 × 35 mm (l × w × h)
- ∴ Increased safety as the AED is readily available
- ∴ Uses the gentle MULTIPULSE BIOWAVE technology



Small enough for  
drone delivery

© Andreas Claesson, Karolinska Institutet, Sweden



## A MASTERPIECE OF SWISS PRECISION

FRED easyport is an extremely small and light pocket defibrillator. This compact lightweight device is the perfect companion for doctors, public service staff and those who always want to have access to an AED. Risk patients have their own rescue device at hand, so that they and their families can react quickly in case of an emergency.

FRED easyport fits in every coat pocket, office drawer, glove compartment or doctor's bag and still fulfils all the demands on a AED.

## FOR ADULTS AND CHILDREN

FRED easyport detects when paediatric electrodes are connected and automatically switches to the paediatric energy level.

Energy settings can be configured for adults as well as children, according to standard values.

## GENTLE ON THE HEART

FRED easyport uses the particularly gentle MULTIPULSE BIOWAVE technology. This technology minimises damage to the myocardium<sup>1</sup> and therefore increases the chances of survival for the patient.



<sup>1</sup> Cansell A. (2000) Wirksamkeit und Sicherheit neuer Impulskurvenformen bei transthorakaler Defibrillation. Biphasische Impulskurvenformen. Notfall & Rettungsmedizin. Springer-Verlag – 3: 458-474. Irnich W, Cansell A. (2005) Principle and Advantage of Chopped Defibrillation Schocks. Biomedizinische Technik – Vol 50 – Suppl. vol. 1 – Part 1: 727 – 728

## CLEAR INSTRUCTIONS

Instructions from loudspeaker and display. FRED easyport is available in different language versions.

## POCKET-SIZE DEFIBRILLATOR


This photo shows its original size (1:1)!

## EASY OPERATION

The 1-2-3 buttons simplify use

## DATA TRANSFER

Mini SD card (option) for ECG and event recording



## LCD DISPLAY

Including ECG waveform display

## DESIGNED TO LAST

Approved for use in ambulances, helicopters and aircrafts

A second battery and a protective bag are included in the scope of delivery.

## FRED EASYPORT TRAINER\*

The FRED easyport is also available as a separate training device. FRED easyport Trainer helps training first responders.

\*no medical device



## MANUAL MODE

The SCHILLER FRED easyport is also available with manual mode. Doctors can switch to manual mode and decide whether or not to deliver a shock.





**Americas**  
**SCHILLER Americas Inc.**  
 Doral, Florida 33172

**North America:**  
 Phone +1 786 845 06 20  
 Fax +1 786 845 06 02  
[sales@schilleramericas.com](mailto:sales@schilleramericas.com)  
[www.schilleramericas.com](http://www.schilleramericas.com)

**Latin America & Caribbean**  
 Phone +1 305 591 11 21  
 Fax +1 786 845 06 02  
[sales@schilleramericas.com](mailto:sales@schilleramericas.com)  
[www.schilleramericas.com](http://www.schilleramericas.com)



**Asia**  
**SCHILLER Asia-Pacific / Malaysia**  
 52200 Kuala Lumpur, Malaysia  
 Phone +603 6272 3033  
[sales@schiller.com.my](mailto:sales@schiller.com.my)  
[www.schiller-asia.com](http://www.schiller-asia.com)



**Austria**  
**SCHILLER Handelsgesellschaft m.b.H.**  
 A-4040 Linz  
 Phone +43 732 70 99 0  
 Fax +43 732 757 000  
[sales@schiller.at](mailto:sales@schiller.at)  
[www.schiller.at](http://www.schiller.at)



**China**  
**Alfred Schiller (Beijing) Medical Equipment Co., Ltd.**  
 100015 Beijing, China  
 Phone +86 010 52007020  
 Fax +86 010 52007020-8016  
[info@schillerchina.com](mailto:info@schillerchina.com)  
[www.schillermedical.cn](http://www.schillermedical.cn)



**Croatia**  
**SCHILLER d.o.o.**  
 10000 Zagreb  
 Phone +385 1 309 66 59  
 Fax +385 1 309 66 60  
[info@schillerzg.hr](mailto:info@schillerzg.hr)  
[www.schiller.hr](http://www.schiller.hr)



**France**  
**SCHILLER Médical**  
 F-67160 Wissembourg  
 Phone +33 3 88 63 36 00  
 Fax +33 3 88 63 36 49  
[info@schiller.fr](mailto:info@schiller.fr)  
[www.schiller.fr](http://www.schiller.fr)



**France (distribution France)**  
**SCHILLER France S.A.S.**  
 F-77608 Bussy St Georges  
 Phone +33 1 64 66 50 00  
 Fax +33 1 64 66 50 10  
[contact@schillerfrance.fr](mailto:contact@schillerfrance.fr)  
[www.schiller-france.com](http://www.schiller-france.com)



**Germany**  
**SCHILLER Medizintechnik GmbH**  
 D-85622 Feldkirchen b. München  
 Phone +49 89 62 99 81 0  
 Fax +49 89 609 50 90  
[info@schillermed.de](mailto:info@schillermed.de)  
[www.schillermed.de](http://www.schillermed.de)



**India**  
**SCHILLER Healthcare India Pvt. Ltd.**  
 Mumbai – 400 059, India  
 Phone +91 22 6152 3333 / 2920 9141  
 Fax +91 22 2920 9142  
[sales@schillerindia.com](mailto:sales@schillerindia.com)  
[www.schillerindia.com](http://www.schillerindia.com)



**Poland**  
**SCHILLER Poland Sp. z o.o.**  
 PL-02-729 Warszawa  
 Phone +48 22 843 20 89 / +48 22 647 35 90  
 Fax +48 22 843 20 89  
[schiller@schiller.pl](mailto:schiller@schiller.pl)  
[www.schiller.pl](http://www.schiller.pl)



**Russia & C.I.S.**  
**AO SCHILLER.RU**  
 119049 Moscow, Russia  
 Phone +7 (495) 970 11 33  
[mail@schiller.ru](mailto:mail@schiller.ru)  
[www.schiller.ru](http://www.schiller.ru)



**Serbia**  
**SCHILLER d.o.o.**  
 11010 Beograd  
 Phone +381 11 39 79 508  
 Fax +381 11 39 79 518  
[info@schiller.rs](mailto:info@schiller.rs)  
[www.schiller.rs](http://www.schiller.rs)



**Slovenia**  
**SCHILLER d.o.o.**  
 2310 Slovenska Bistrica  
 Phone +386 2 843 00 56  
 Fax +386 2 843 00 57  
[info@schiller.si](mailto:info@schiller.si)  
[www.schiller.si](http://www.schiller.si)



**Spain**  
**SCHILLER ESPAÑA, S.A.**  
 ES-28232 Las Rozas/Madrid  
 Phone +34 91 713 01 76  
 Fax +34 91 355 79 33  
[schiller@schiller.es](mailto:schiller@schiller.es)  
[www.schiller.es](http://www.schiller.es)



**Switzerland**  
**SCHILLER-Reomed AG**  
 CH-8912 Obfelden  
 Phone +41 44 744 30 00  
 Fax +41 44 740 37 10  
[sales@schiller-reomed.ch](mailto:sales@schiller-reomed.ch)  
[www.schiller-reomed.ch](http://www.schiller-reomed.ch)



**Turkey**  
**SCHILLER TÜRKIYE**  
 Okmeydanı-Sisli – Istanbul  
 Phone +90 212 210 8681 (pbx)  
 Fax +90 212 210 8684  
[info@schiller.com.tr](mailto:info@schiller.com.tr)  
[www.schiller-turkiye.com](http://www.schiller-turkiye.com)



**United Kingdom**  
**SCHILLER UK Ltd.**  
 Anstruther, Fife KY10 3H  
 Phone +44 16 190 72 565  
[sales@schilleruk.com](mailto:sales@schilleruk.com)  
[www.schilleruk.com](http://www.schilleruk.com)

Device availability in your market is subject to regulatory approval.

