




## FRED easy Life

The Automated External Defibrillator that gives everyone the confidence to save lives.



**SCHILLER**

The Art of Diagnostics

# FRED easy Life

**The Automated External Defibrillator that gives everyone the confidence to save lives.**

Simple, clear and detailed voice and displayed instructions are complemented by pictures on the front panel to give visual reminders to the rescuer.

## Guided from the beginning to the end

- Pictures on the front panel remind you of the main resuscitation steps.
- Voice and display instructions guide you in all environmental conditions (noisy, dark).
- A metronome gives you the correct rhythm to perform good chest compressions.

# Anybody can save



1 Call emergency services



4 Do not touch the patient  
Follow instructions



## Anybody can save everybody...

### Suitable for all users, available in:

- Fully-automatic mode: attach electrodes to the patient, switch the unit on and let the device analyse and deliver the shock if necessary.
- Semi-automatic mode: attach electrodes to the patient, switch the unit on and follow the instructions, press the shock button if recommended by the device.

### Designed for all patients:

Each patient is different. Thanks to automatic impedance measurement, the device adapts the treatment to each pa-

tient in order to reach the same high level of efficacy.

FRED easy Life is adapted to adults and pediatrics: electrodes are automatically recognised and the energy level is adjusted accordingly.

## ...everywhere

Around the world, public places, schools, companies, airports, stadiums and shopping centres equip themselves with Automated External Defibrillators (AED), to guarantee security and first aid for everyone.

With its new concept of networked AEDs, SCHILLER advises a unique solution for a better reactivity in case of cardiac arrest and a simplified management of the defibrillator pool.

FRED® easy Life Online is an AED which is permanently under remote control (presence, functioning, configuration view and update...).

Save everybody ...everywhere.



2 Switch unit on



3 Place electrodes on chest



6 Perform chest compressions

5 Deliver shock if necessary

**Available versions:**

FRED easy Life fully automatic



FRED easy Life semi automatic




FRED easy Life Online fully automatic



FRED easy Life Online semi automatic



**Example of a typical installation controlled via the internet:**







**Asia**  
**SCHILLER Asia-Pacific / Malaysia**  
52200 Kuala Lumpur, Malaysia  
Phone +603 6272 3033  
Fax +603 6272 2030  
sales@schiller.com.my  
www.schiller-asia.com



**Austria**  
**SCHILLER Handelsgesellschaft m.b.H.**  
A-4040 Linz  
Phone +43 732 709 90  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at



**China**  
**Alfred Schiller Medical Equipment Co. Ltd.**  
100015 Beijing, China  
Phone +86-010-52007020  
diamond@schillermedical.com  
www.schiller.cn



**France**  
**SCHILLER Médical S.A.S.**  
F-67162 Wissembourg/Cedex  
Phone +33 3 88 63 36 00  
Fax +33 3 88 94 12 82  
info@schiller.fr  
www.schiller-medical.com



**France (distribution France)**  
**SCHILLER France S.A.S.**  
F-77600 Bussy St Georges  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
infoschiller@schiller-france.fr  
www.schiller-france.com



**Germany**  
**SCHILLER Medizintechnik GmbH**  
D-85622 Feldkirchen b. München  
Phone +49 89 62 99 81-0  
Fax +49 89 62 99 81-54  
info@schillermed.de  
www.schillermed.de



**Hungary**  
**SCHILLER Diamed Ltd.**  
H-1141 Budapest  
Phone +36 (1) 383-4780 / 460-9491  
Fax +36 (1) 383-4778  
sales@schiller.at  
www.schiller-hungary.hu



**India**  
**SCHILLER Healthcare India Pvt. Ltd.**  
Mumbai - 400 001, India  
Phone +91 22 6152 3333/ 2920 9141  
Fax +91 22 2920 9142  
sales@schillerindia.com  
www.schillerindia.com



**Japan**  
**SCHILLER Japan, Ltd.**  
Hiroshima 734-8551  
Phone +81 82 250 2055  
Fax +81 82 253 1713  
koji.maekawa@schiller.jp  
www.schiller.jp



**Croatia**  
**Schiller medicinski instrumenti d.o.o.**  
10000 Zagreb  
Phone +385 1 309 66 59  
Fax +385 1 309 66 60  
info@schillerzg.hr  
www.schiller.ch



**Latin America**  
**SCHILLER Latin America, Inc.**  
Doral, Florida 33172  
Phone +1 954 673 0358  
Fax +1 786 845 06 02  
info@schillerla.ch  
www.schillerla.ch



**Poland**  
**SCHILLER Poland Sp. z o.o.**  
PL-02-729 Warszawa  
Phone +48 22 8432089  
Fax +48 22 8432089  
schiller@schiller.pl  
www.schiller.pl



**Russia & C.I.S.**  
**SCHILLER AG Rep. office**  
125124 Moscow, Russia  
Phone +7 (495) 970 11 33  
Fax +7 (495) 970 11 33  
mail@schiller-ag.com  
www.schiller-cis.com



**Serbia**  
**SCHILLER d.o.o.**  
11010 Beograd  
Phone +381 11 39 79 508  
Fax +381 11 39 79 518  
info@schiller.rs  
www.schiller.rs



**Slovenia**  
**SCHILLER d.o.o.**  
2310 Slovenska Bistrica  
Phone +386 2 843 00 56  
Fax +386 2 843 00 57  
info@schiller.si  
www.schiller.si



**Spain**  
**SCHILLER ESPAÑA, S.A.**  
E-28230-Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller.es



**Switzerland**  
**SCHILLER-Reomed AG**  
CH-8953 Dielikon  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch



**Turkey**  
**SCHILLER TÜRKIYE**  
Okmeydanı-Sisli - İstanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
sales@schiller-turkiye.com  
www.schiller-turkiye.com



**USA**  
**SCHILLER America Inc.**  
Doral, Florida 33172  
Phone +1 786 845 0620  
Fax +1 786 845 06 02  
sales@schilleramerica.com  
www.schilleramerica.com

08/11

Part no.: 2.500610



# SCHILLER

The Art of Diagnostics