



# CARDIOVIT CS-200 Office

Wireless ECG acquisition




# CARDIOVIT CS-200 OFFICE

The CARDIOVIT CS-200 Office offers cutting-edge diagnostic capability for 12-lead resting and stress wireless ECG acquisition. The choice is yours: Whether only the software for any laptop or as a package integrated in a medical PC, it's the perfect solution for GPs, clinics and hospitals.



MS-12 BLUE



**SWITCH BETWEEN TWO CUSTOMISABLE  
SCREEN LAYOUTS AT THE TOUCH OF A  
BUTTON**

## ECG ACQUISITION WITH MAXIMUM FREEDOM

With the introduction of the MS-12 blue, SCHILLER provides a new solution that greatly simplifies ECG acquisition.

- ❖ Allows to freely operate in any environment without worrying about cable length.
- ❖ Weighing only 115g and fitting in one hand, the MS-12 blue sits on the bed or on the patient body during ECG acquisition.
- ❖ Motion artefact is minimized without an analogue cable connection and an extremely lightweight patient cable.
- ❖ One battery charge guarantees continuous operation for more than 36 hours and more than a week with the power-saving mode.

## ADVANCED CLINICAL TOOLS

The CS-200 Office can be supplemented with high-end software features:

- ❖ ETM interpretation
- ❖ ETM sport
- ❖ 12-lead exercise ECG
- ❖ 3D vector ECG
- ❖ Arrhythmia detection and classification
- ❖ Risk factors

## USER-DEFINED WORKFLOW (WITH PRINTOUTS OR PAPERLESS)

The CARDIOVIT CS-200 Office fulfills the needs of every user thanks to a customizable user interface, stage and report printout with laser printer, and seamless connectivity to EMR, PACS, HIS or SCHILLER SEMA data management system.

## SOLUTIONS FOR ANY ENVIRONMENT

The CS-200 Office can be easily installed on your laptop or PC, as well as purchased as a ready-to-use system including a robust trolley and a 21.5" touch panel PC with pre-installed software (see picture).

Note:  
CARDIOVIT CS-200 Office is also available with MS-12 USB (USB connection).

High-resolution and easy-to-clean 21.5" touch screen.

Intuitive and clear user interface suitable for touch screen.

Soft handle positioned at an adjustable height.

Support for PC power supply and cables.

Transparent bin to store and transport medical items.

User-defined screen layouts and workflow.

MS-12 blue: Small but powerful wireless ECG acquisition module with on-screen signal view and hookup advisor.

Easy-glide wheels to smoothly maneuver the trolley.





**Americas**  
SCHILLER Americas Inc.  
Doral, Florida 33172

**North America:**  
Phone +1 786 845 06 20  
Fax +1 786 845 06 02  
sales@schilleramerica.com  
www.schilleramerica.com

**Latin America & Caribbean**  
Phone +1 305 591 11 21  
Fax +1 786 845 06 02  
sales@schillerla.ch  
www.schillerla.ch



**Asia**  
SCHILLER Asia-Pacific / Malaysia  
52200 Kuala Lumpur, Malaysia  
Phone +603 6272 3033  
Fax +603 6272 2030  
sales@schiller.com.my  
www.schiller-asia.com



**Austria**  
SCHILLER Handelsgesellschaft m.b.H.  
A-4040 Linz  
Phone +43 732 70 99  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at



**China**  
Alfred Schiller (Beijing) Medical Equipment  
Co., Ltd.  
100015 Beijing, China  
Phone +86 010 52007020  
Fax +86 010 52007020-8016  
info@schillerchina.com  
www.schillermedical.cn



**Croatia**  
SCHILLER d.o.o.  
10000 Zagreb  
Phone +385 1 309 66 59  
Fax +385 1 309 66 60  
info@schillerzg.hr  
www.schiller.hr



**France**  
SCHILLER Médical S.A.S.  
F-67162 Wissembourg  
Phone +33 3 88 63 36 00  
Fax +33 3 88 63 36 49  
info@schiller.fr  
www.schiller-medical.com



**France (distribution France)**  
SCHILLER France S.A.S.  
F-77608 Bussy St Georges  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
contact@schillerfrance.fr  
www.schiller-france.com



**Germany**  
SCHILLER Medizintechnik GmbH  
D-85622 Feldkirchen b. München  
Phone +49 89 62 99 81 0  
Fax +49 89 609 50 90  
info@schillermed.de  
www.schillermed.de



**Hungary**  
SCHILLER Diamed Ltd.  
H-1141 Budapest  
Phone +36 1 383 4780 / +36 1 460 9491  
Fax +36 1 383 4778  
info@schillerhungary.hu  
www.schillerhungary.hu



**India**  
SCHILLER Healthcare India Pvt. Ltd.  
Mumbai – 400 059, India  
Phone +91 22 6152 3333 / 2920 9141  
Fax +91 22 2920 9142  
sales@schillerindia.com  
www.schillerindia.com



**Poland**  
SCHILLER Poland Sp. z o.o.  
PL-02-729 Warszawa  
Phone +48 22 843 20 89 / +48 22 647 35 90  
Fax +48 22 843 20 89  
schiller@schiller.pl  
www.schiller.pl



**Russia & C.I.S.**  
AO SCHILLER.RU  
125124 Moscow, Russia  
Phone +7 (495) 970 11 33  
Fax +7 (495) 956 29 10  
mail@schiller.ru  
www.schiller.ru  
www.schiller-cis.com



**Serbia**  
SCHILLER d.o.o.  
11010 Beograd  
Phone +381 11 39 79 508  
Fax +381 11 39 79 518  
info@schiller.rs  
www.schiller.rs



**Slovenia**  
SCHILLER d.o.o.  
2310 Slovenska Bistrica  
Phone +386 2 843 00 56  
Fax +386 2 843 00 57  
info@schiller.si  
www.schiller.si



**Spain**  
SCHILLER ESPAÑA, S.A.  
ES-28232 Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller.es



**Switzerland**  
SCHILLER-Reomed AG  
CH-8953 Dietikon  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch



**Turkey**  
SCHILLER TÜRKİYE  
Okmeydanı-Sisli – Istanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
info@schiller.com.tr  
www.schiller-turkiye.com



**United Kingdom**  
SCHILLER UK  
Bellshill, ML4 3PR  
Phone +44 1698 744 505  
Fax +44 1698 744 474  
sales@schilleruk.com  
www.schilleruk.com

