

SCHILLER

Home

Patient

DIAGNOSTIC STATION DS20

Modern and practical design unites
the most demanded functions

BP
Hg

Start BP

100

120

120 / 80

160

200

(92)

Mode: Manual

07:55

HR
bpm

80

SpO2
% Normal

SpCO
%

2

RR
rpm

15

TEMP
°C


37.0

07:59

WEIGHT
kg

72.6

07:55





SCHILLER

The Art of Diagnostics

DIAGNOSTIC STATION DS20

SCHILLER's DIAGNOSTIC STATION DS20 simplifies your daily work:

- ✦ Main vital signs and physical assessment tools in one device
- ✦ Ideal for arrhythmia and atrial fibrillation screening: one ECG channel for display, storage and printing, with only 3 electrodes and without undressing the patient
- ✦ Large interactive touch screen and self-explanatory user interface for maximum user-friendliness



FAST AND EASY TO USE

The DS20 immediately detects connected sensors and automatically displays the corresponding value. This fast operation allows for easy addition or removal of parameters.

Thanks to its intuitive user guidance, SCHILLER's DS20 is self-explanatory and very little training is needed. The large interactive touch screen supports the highest user-friendliness.

ONE TOUCH TO THE MEASUREMENT

Simply touch the large, high-resolution colour display and you are in business recording, and selecting the highest quality ECGs.

- ✦ 1. Perform
- ✦ 2. Review
- ✦ 3. Store or export wherever you want

CONNECTED

The DIAGNOSTIC STATION DS20 is a networked device. Seamless connectivity to EMR, PACS, HIS or SCHILLER's SEMA3 Cardiology Information System is possible and bidirectional communication allows for easy data access, while Wi-Fi with strong security enables direct and fast transmission.

EXPANDABLE

Easily add new functions and other devices or future technological developments. The DS20 is ready for the most common functions and will satisfy new requirements as they arise in your institution.

See and touch 18.5-inch responsive and high-resolution multi-touch screen.

The DS20 offers optional 12-lead resting ECG with a sophisticated interpretation, as well as resting rhythm

Full Spirometry testing (FVC, SVC, Pre/Post)

Pulse oximetry: choice of Masimo or Nellcor technology. Masimo Rainbow SpCO parameter optionally available.

Non-invasive temporal artery thermometry for fast, accurate temperature with a gentle forehead scan

"Pulse Wave Analysis" (PWA), Analysis of pulse wave velocity and blood pressure measurement in one





Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 99
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment
Co., Ltd.
100015 Beijing, China
Phone +86 010 52007020
Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical S.A.S.
F-67162 Wissembourg
Phone +33 3 88 63 36 00
Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 1 383 4780 / +36 1 460 9491
Fax +36 1 383 4778
info@schillerhungary.hu
www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 059, India
Phone +91 22 6152 3333 / 2920 9141
Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru
www.schiller-cis.com



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
Okmeydanı-Sisli - Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
Bellshill, ML4 3PR
Phone +44 1698 744 505
Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com

Article no.: 2.500692
2017/06_EN



SCHILLER
The Art of Diagnostics