

# SCOUTTUBE

## THE NEW STANDARD IN SPIROMETRY

ScoutTubes ensure a safe, accurate and patient friendly spirometry test with the SpiroScout, the SpiroScout SP and the SpiroScout SP plus devices.

The main advantages of the ScoutTubes are:

- ✦ Always accurate
- ✦ Resistance free
- ✦ Calibration free
- ✦ Disposable, therefore no cleaning
- ✦ Maintenance free
- ✦ No errors caused by miscalibration

- ✦ Ensures a tight seal around the lips
- ✦ Designed for all mouth shapes and sizes
- ✦ Suits all ages from children to adults, at all levels of fitness



### POSITIONING OF TEETH

- ✦ Keeps ScoutTube securely in place

### TIGHT LIP SEAL

- ✦ Prevents leaks

### MOUTHPIECE

- ✦ Shaped oval for optimal lip seal and teeth positioning
- ✦ Conical to suit different mouth sizes

### FILTER MEMBRANE

- ✦ Bacteria filter membrane built into each ScoutTube
- ✦ Prevents cross-contamination





# Benefits of ScoutTube

## ECONOMICAL

- ∴ No running costs
- ∴ No costs for calibration equipment
- ∴ No costs or downtime for maintenance

## SIMPLE APPLICATION

- ∴ Just insert the ScoutTube, ready for measurement

## MAINTENANCE-FREE

- ∴ No maintenance is needed with ScoutTube or SpiroScout
- ∴ No calibration or checking procedures are required

## HYGIENIC

- ∴ One of the safest respiratory flow measurement systems available
- ∴ Can be fitted with the GANSHORN PFT filter

## ENVIRONMENTALLY SAFE

- ∴ Lightweight PE, environmentally friendly, disposable and bio-degradable

## DISPOSABLE AND SUSTAINABLE

To ensure a sustainable product, SCHILLER compensates the CO<sub>2</sub> emitted during production, packaging, shipment and disposal of the ScoutTubes and the PFT filters.



ScoutTube is compatible with the following devices:

- ∴ SpiroScout
- ∴ SpiroScout SP
- ∴ SpiroScout SP plus

Subject to changes.  
2027/05, EN



Article no.: 2.500743



**SCHILLER**  
The Art of Diagnostics